


Introducing the Anaesthetic CADDY



Does your operating room look like this?



Effortlessly transform it into this!



“Optimise your operating room efficiency with the Anaesthetic CADDY”



Efficient Workflow

Holds anaesthetic circuit tubing, in-line processing modules, sucker tubing, and monitoring cables in designated slots.



Easy Cleaning

Completely reusable & compatible with standard hospital cleaning agents.



Secure Attachment

Firmly attaches to the operating table with Velcro for stability. Strong and durable.




Eco-Friendly

Made from 100% recyclable polycarbonate. Latex Free




Anaesthetic Caddy & Arterial Line Holder



Anaesthetic Caddy

Designed for anaesthetic use to hold airway tubing, suction tubing, monitoring modules, spare cabling and much more.


Width: 300 mm
SKU: AC-PC-V3-001



Anaesthetic Caddy - Beach Chair

Designed for anaesthetic use when the patient is in the beach chair position. This Caddy can hold airway tubing, suction tubing, monitoring modules, spare cabling and much more.

Width: 200 mm
SKU: ACBC-PC-V1-001



Arterial Line Holder

Designed to attach to the Anaesthetic Caddy. The arterial line holder is engineered to hold the arterial line and reposition airway tubing to an elevated position.

Width: 300 mm
SKU: ACAL-PC-V1-001

1300 144 734

sales@lateralmedical.com
www.lateralmedical.com